

Wound Care < 40 YO & Never Bent

IO Takes 1st 10 of 2nd O₂

Hyperbaric Chamber Profile


7 Min descent: IO SE = 56, Repet = F

9 Min descent: IO SE = 56, Repet = F

11 Min descent: IO SE = 55, Repet = F

13 Min descent: IO SE = 55, Repet = F

15 Min descent: IO SE = 54, Repet = F



Wound Care > 39 or Bent

IO Takes 1st 15 of 2nd O₂

Hyperbaric Chamber Profile


7 Min descent: IO SE = 62, Repet = E

9 Min descent: IO SE = 62, Repet = E

11 Min descent: IO SE = 61, Repet = E

13 Min descent: IO SE = 61, Repet = E

15 Min descent: IO SE = 60, Repet = E



Wound Care < 40 YO & Never Bent

IO Takes Last 10 of 2nd O₂

Hyperbaric Chamber Profile


7 Min descent: IO SE = 56, Repet = F

9 Min descent: IO SE = 56, Repet = F

11 Min descent: IO SE = 55, Repet = F

13 Min descent: IO SE = 55, Repet = F

15 Min descent: IO SE = 54, Repet = F



Wound Care > 39 or Bent

IO Takes 2nd 15 of 2nd O₂

Hyperbaric Chamber Profile


7 Min descent: IO SE = 62, Repet = E

9 Min descent: IO SE = 62, Repet = E

11 Min descent: IO SE = 61, Repet = E

13 Min descent: IO SE = 61, Repet = E

15 Min descent: IO SE = 60, Repet = E



Wound Care: 2 IOs Split

Each IO Takes 2x10 Min

Patients O2 shown in dark orange

